

Product specifications

Name	Anti-Influenza B 9906 SPTN-5
Specificity	Antibody recognizes <i>Influenza</i> B nucleoprotein (np)
Description	Monoclonal mouse antibody, cultured <i>in vitro</i> under conditions free from animal-derived components.
Product code	100437
Product buffer solution	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative
Shelf life and storage	Unspecified, storage at 2–8 °C
Subclass	IgG _{2b}
Analyte description	Influenza B virus is the only species in a genus in the virus family Orthomyxoviridae. Influenza B viruses are only known to infect humans and seals. This limited host range is apparently responsible for the lack of Influenza B caused influenza pandemics in contrast with those caused by the morphologically similar Influenza A. In addition, Influenza B mutates at a rate 2-3 times lower than type A. However, influenza B mutates enough that lasting immunity is not possible.

Parameters tested on each lot

Product appearance	Liquid, may turn slightly opaque during storage
Product concentration	5.0 mg/ml (+/- 10 %)
Immunoreactivity	80–120 % compared to the reference sample in an FIA test
IEF Profile	6.8–7.4
Purity	≥ 95 %

Kinetic parameters

Association rate constant	7.6×10^4 1/Ms
Dissociation rate constant	8.0×10^{-4} 1/s
Affinity constant	$K_A = 9.5 \times 10^7$ 1/M; $K_D = 1.1 \times 10^{-8}$ M (=11.0 nM)
Determination method	SPR analysis (ProteOn XPR36)
Determination antigen	Recombinant Influenza B virus nucleoprotein from strain B/Singapore/222/79 (in-house antigen, UniProtKB P04666)

Legal disclaimer

Medix Biochemica products meet their specifications if transported, stored and used according to the instructions. MedixMAB and MedixAntigens are trademarks of Medix Biochemica and may not be used or reproduced without Medix Biochemica's written permission.



Cross-reactivities	Recognizes Influenza B/Yamagata clade 3, Influenza B/Victoria clade 1A, clade 1A ($\Delta 2$), and clade 1A ($\Delta 3$). Does not recognize Influenza A (H3N2) group 3C.2a1b and group 3C.3a, Influenza A (H1N1)pdm09 group 6B.1A1 and group 6B.1A5.																
Epitope	Not determined (N/D)																
Pair recommendations	As nucleoprotein is abundant in the virus, all Anti-Influenza B antibodies (9901, 9905, 9906, 9908 and 9909) pair with themselves and each other.																
Platforms tested	FIA																
Antigens tested	Recombinant Influenza B NP antigen, Medix Biochemica 610051																
Product stability	<table><thead><tr><th>TEMPERATURE, TIME</th><th>RESULT</th></tr></thead><tbody><tr><td>-70 °C, 21 days</td><td>OK</td></tr><tr><td>-20 °C, 21 days</td><td>OK</td></tr><tr><td>+4 °C, 21 days</td><td>OK</td></tr><tr><td>+35 °C, 7 days</td><td>OK</td></tr><tr><td>+35 °C, 21 days</td><td>Some charge alterations</td></tr><tr><td>+45 °C, 3 days</td><td>OK</td></tr><tr><td>+45 °C, 7 days</td><td>Some charge alterations</td></tr></tbody></table> <p>Stability testing is performed in the product buffer to see whether different temperatures affect the antigen binding, charge or composition of the antibody. Please note that the shelf life given on the first page is based on real time stability testing at 2–8 °C in the product buffer.</p>	TEMPERATURE, TIME	RESULT	-70 °C, 21 days	OK	-20 °C, 21 days	OK	+4 °C, 21 days	OK	+35 °C, 7 days	OK	+35 °C, 21 days	Some charge alterations	+45 °C, 3 days	OK	+45 °C, 7 days	Some charge alterations
TEMPERATURE, TIME	RESULT																
-70 °C, 21 days	OK																
-20 °C, 21 days	OK																
+4 °C, 21 days	OK																
+35 °C, 7 days	OK																
+35 °C, 21 days	Some charge alterations																
+45 °C, 3 days	OK																
+45 °C, 7 days	Some charge alterations																
Miscellaneous	-																
References	-																

