Medix Biochemica

Anti-h NT-proBNP 1310 SPTN-5

Product overview					
Catalog number	100718				
Specificity	Antibody recognizes human N-terminal pro-B-type natriuretic peptide				
Description	Monoclonal mouse antibody, cultured <i>in vitro</i> under conditions free from animal-derived components				
Product buffer solution	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN $_3$ as a preservative				
Shelf life and storage	Unspecified, storage at 2–8 °C				
Subclass	IgG1				
Analyte description	The plasma concentrations of both the N-terminal fragment of Pro-B-type natriuretic peptide (NT-proBNP) and prohormone of brain natriuretic peptide (proBNP) are significantly increased in patients with asymptomatic and symptomatic left ventricular dysfunction.				
Parameters tested on each lot					
Product appearance	Liquid, may turn slightly opaque during storage				
Product concentration	5.0 mg/ml (+/-10 %)				
Immunoreactivity	80–120 % compared to the reference sample in an FIA test				
IEF Profile	6.6–6.9				
Purity	≥ 95 %				

Kinetic parameters	
Association rate constant	2.0 x 10 ⁶ 1/Ms
Dissociation rate constant	3.4 x 10 ⁻⁴ 1/s
Affinity constant	$K_A = 5.8 \times 10^9 \text{ 1/M}; K_D = 1.9 \times 10^{-10} \text{ M} (= 0.19 \text{ nM})$
Determination method	BLI (Octet RED96e)
Determination antigen	Recombinant NT-proBNP, Medix Biochemica (Cat 610090)



Oy Medix Biochemica Ab Klovinpellontie 3, Fl-02180 Espoo, Finland medix@medixbiochemica.com www.medixbiochemica.com VAT reg. no. FI14631532

Legal disclaimer

Medix Biochemica products meet their specifications if transported, stored, and used according to the instructions. Medix Biochemica's products may not be used or reproduced without Medix Biochemica's written permission.

Medix Biochemica

Cross-reactivities Does not recognize canine NT-proBNP (HyTest, 8CNT9).

Epitope

Amino acids 57-76 (ATEGIRGHRKMVLYTLRAPR)

Pair recommendation

ns			DETECTION									
			1306/ RC1306	1307/ RC1307	1308	1309	1310	1311	1312	HM071	HM148	HM210
		1306 / RC1306	-	+	-	+	+	+	+	-	-	-
		1307/ RC1307	+	-	+	+	+	+	+	-	+	-
		1308	-	+	-	+	+	+	+	-	-	-
	RE	1309	+	+	+	-	-	-	-	-	-	-
	CAPTURE	1310	+	+	+	-	-	-	-	+	-	+
	ç	1311	+	+	+	-	-	-	-	-	-	-
		1312	+	+	+	-	-	-	-	-	-	-
		HM071	-	-	-	+	-	+	-	-	-	-
		HM146	-	-	-	-	+	+	+	-	+	-

Following pairs are especially recommended for the below mentioned assays: CLIA: 1306 (capture) - 1309 (detection), 1306 - 1311, 1306 - 1312, 1308 - 1309, 1308 - 1310, 1308 - 1311 and 1308 - 1312 LF: 1306 (capture) - 1309 (detection), 1306 - 1310 FIA: 1306 (capture) - 1310 (detection), 1306 - 1312, 1309 - 1307, 1310 - HM210, HM071 - 1309, 1310 - HM071, 1307 - 1311

Please note that pair recommendations are based on results obtained by our laboratory. Equally good results may be obtained using other pairs and therefore these recommendations are only indicative.

Platforms tested FIA, LF, CLIA

Antigens tested Recombinant NT-proBNP antigen, Medix Biochemica 610090 and recombinant proBNP antigen Medix Biochemica 710017, 710042 and 710043.

Product stability	TEMPERATURE, TIME	RESULT
	-70 °C, 21 days	OK
	-20 °C, 21 days	OK
	+4 °C, 21 days	OK
	+35 °C, 21 days	OK
	+45 °C, 3 days	OK
	+45 °C, 7 days	Reduced homogeneity

Stability testing is performed in the product buffer to see whether different temperatures affect the antigen binding, charge or composition of the antibody. Please note that the shelf life given on the first page is based on real time stability testing at 2–8 °C in the product buffer.

```
Miscellaneous
```

References



Oy Medix Biochemica Ab Klovinpellontie 3, FI-02180 Espoo, Finland

medix@medixbiochemica.com www.medixbiochemica.com VAT reg. no. FI14631532

Legal disclaimer

Medix Biochemica products meet their specifications if transported, stored, and used according to the instructions. Medix Biochemica's products may not be used or reproduced without Medix Biochemica's written permission.