# Medix Biochemica

## Anti-h BNP 11908 SPTN-5

| Product overview              |  |  |  |  |
|-------------------------------|--|--|--|--|
| Catalog number                | 100973   |  |  |  |
| Specificity                   | Antibody recognizes human B-type natriuretic peptide   |  |  |  |
| Description                   | Monoclonal mouse antibody, cultured <i>in vitro</i> under conditions free from animal-derived components.  |  |  |  |
| Product buffer solution       | 50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN $_3$ as a preservative   |  |  |  |
| Shelf life and storage        | Unspecified, storage at 2–8 °C   |  |  |  |
| Subclass                      | lgG <sub>1</sub>   |  |  |  |
| Analyte description           | B-type natriuretic peptide (BNP) is a cardiac hormone released by the<br>heart in response to ventricular myocardium wall stress. In patients with<br>heart failure, BNP levels are elevated and assessed as important<br>measures of cardiac function and diagnosis of heart failure. |  |  |  |
| Parameters tested on each lot |  |  |  |  |
| Product appearance            | Liquid, may turn slightly opaque during storage  |  |  |  |
| Product concentration         | 5.0 mg/ml (+/- 10%)  |  |  |  |
| Immunoreactivity              | 80–120% compared to the reference sample in an FIA test  |  |  |  |
| IEF Profile                   | 6.1–6.5  |  |  |  |
| Purity                        | ≥ 95 %   |  |  |  |
| Kinetic parameters            |  |  |  |  |
| Association rate constant     | 1.1 x 10 <sup>6</sup> 1/Ms   |  |  |  |

**Dissociation rate constant** Does not dissociate under conditions used.

Affinity constantN/A (Not applicable)

- Determination method BLI (Octet RED96e)
- Determination antigen Synthetic BNP, Proteogenix



**Oy Medix Biochemica Ab** Klovinpellontie 3, Fl-02180 Espoo, Finland medix@medixbiochemica.com www.medixbiochemica.com VAT reg. no. FI14631532

#### Legal disclaimer

Medix Biochemica products meet their specifications if transported, stored, and used according to the instructions. Medix Biochemica's products may not be used or reproduced without Medix Biochemica's written permission.

# Medix Biochemica

### **Cross-reactivities**

Does not recognize ANP or CNP.

Amino acid region 10-21.

Epitope

Pair recommendations

|         |       | DETECTION |       |       |
|---------|-------|-----------|-------|-------|
|         |       | 11904     | 11906 | 11908 |
| CAPTURE | 11904 | -         | +     | +     |
|         | 11906 | +         | -     | -     |
|         | 11908 | +         | -     | -     |

Please note that pair recommendations are based on results obtained by our laboratory. Equally good results may be obtained using other pairs and therefore these recommendations are only indicative.

| Platforms tested  | FIA  |  |  |
|-------------------|--|--|--|
| Antigens tested   | Synthetic Brain Natriuretic Peptide-32 (BNP), Medix Biochemica 129-10 and synthetic BNP antigen, Medix Biochemica LA788. |  |  |
| Product stability | temperatures affect the antigen bind   | life given on the first page is based on |  |
| Miscellaneous     | -  |  |  |
| References        | -  |  |  |



**Oy Medix Biochemica Ab** Klovinpellontie 3, Fl-02180 Espoo, Finland

medix@medixbiochemica.com www.medixbiochemica.com VAT reg. no. FI14631532

#### Legal disclaimer

Medix Biochemica products meet their specifications if transported, stored, and used according to the instructions. Medix Biochemica's products may not be used or reproduced without Medix Biochemica's written permission.