

2023-12-11



Anti-h PIGF NA000302 SPTN-1

Product overview

Catalog number C-10-0021

Specificity Antibody recognizes human PIGF (Placental growth factor)

Description Recombinant chimeric antibody containing variable region from mouse

IgG and constant region from human IgG1, cultured in vitro under

conditions free from animal-derived components.

Product buffer solution 50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.1% ProClinTM 300 as a preservative

Shelf life and storage Unspecified, storage at 2–8 °C

Subclass Human IgG₁

Analyte description PIGF is a member of vascular endothelial growth factor (VEGF) family. Its

molecular structure is glycoprotein homodimer, which has four isoforms. PIGF stimulate endothelial cell proliferation and migration and it binds to the receptor FLT1/VEGFR-1. PIGF is always used in combination with sFlt-1 to determine the sFlt-1/PIGF ratio. The sFlt-1/PIGF ratio is intended for use as an aid in the diagnosis of preeclampsia in conjunction with other diagnostic

and clinical information.

Parameters tested on each lot

Product appearance Liquid, may turn slightly opaque during storage

Product concentration 1.0 mg/ml (+/- 10%)

Immunoreactivity 80–120% compared to the reference sample

IEF Profile

Purity ≥ 95 %

Kinetic parameters

Association rate constant Not Determined (N/D)

Dissociation rate constant N/D

Affinity constant N/D

Determination method -

Determination antigen -



上海市闵行区浦江绿洲环路10号6幢11层

电话: 021 68119180, 68119181, 68119105

邮箱: medixchina@medixbiochemica.com

网址: www.cnmedixbiochemica.com



Product specification ANTIBODY

2023-12-11

Cross-reactivities

Epitope N/D

Pair recommendations

		DETECTION			
		NA000301	NA000302	NA000303	NA000305
CAPTURE	NA000301	-	-	+	-
	NA000302	-	-	+	-
	NA000303	+	+	-	+
	NA000305	-	-	+	-

Please note that pair recommendations are based on results obtained by our laboratory. Equally good results may be obtained using other pairs and therefore these recommendations are only indicative.

Platforms tested ELISA, CLIA

Antigens tested -

Product stability TEMPERATURE, TIME RESULT

-70 °C, 21 days OK
-20 °C, 21 days OK
+4 °C, 21 days OK
+35 °C, 21 days OK
+45 °C, 7 days OK

Stability testing is performed in the product buffer to see whether different temperatures affect the antigen binding, charge or composition of the antibody. Please note that the shelf life given on the first page is based on real time stability testing at 2–8 °C in the product buffer.

电话: 021 68119180, 68119181, 68119105 邮箱: medixchina@medixbiochemica.com

网址: www.cnmedixbiochemica.com

Miscellaneous -

References -

