# Medix Biochemica

## Anti-h BNP 11906 SPTN-5

Product overview				
Catalog number	100971			
Specificity	Antibody recognizes human B-type natriuretic peptide			
Description	Monoclonal mouse antibody, cultured <i>in vitro</i> under conditions free from animal-derived components.			
Product buffer solution	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN $_3$ as a preservative			
Shelf life and storage	Unspecified, storage at 2–8 °C			
Subclass	lgG1			
Analyte description	B-type natriuretic peptide (BNP) is a cardiac hormone released by the heart in response to ventricular myocardium wall stress. In patients with heart failure, BNP levels are elevated and assessed as important measures of cardiac function and diagnosis of heart failure.			
Parameters tested on each lot				
Product appearance	Liquid, may turn slightly opaque during storage			
Product concentration				
	5.0 mg/ml (+/- 10%)			
Immunoreactivity	5.0 mg/ml (+/- 10%) 80–120% compared to the reference sample in an FIA test			
Immunoreactivity IEF Profile				
-	80–120% compared to the reference sample in an FIA test			
IEF Profile	80–120% compared to the reference sample in an FIA test 7.2–8.0			
IEF Profile Purity	80–120% compared to the reference sample in an FIA test 7.2–8.0			

**Affinity constant**  $K_A = 9.9 \times 10^9 \text{ 1/M}; K_D = 1.0 \times 10^{-10} \text{ M} (= 0.10 \text{ nM})$ 

- Determination method BLI (Octet RED96e)
- Determination antigen Synthetic BNP, Proteogenix



**Oy Medix Biochemica Ab** Klovinpellontie 3, Fl-02180 Espoo, Finland medix@medixbiochemica.com www.medixbiochemica.com VAT reg. no. FI14631532

#### Legal disclaimer

Medix Biochemica products meet their specifications if transported, stored, and used according to the instructions. Medix Biochemica's products may not be used or reproduced without Medix Biochemica's written permission.

# Medix Biochemica

### **Cross-reactivities**

Does not recognize ANP or CNP.

Epitope

Pair recommendations

N/D (N-terminal or mid-region based on pairing information)

		DETECTION		
		11904	11906	11908
APTURE	11904	-	+	+
	11906	+	-	-
CA	11908	+	-	-

Please note that pair recommendations are based on results obtained by our laboratory. Equally good results may be obtained using other pairs and therefore these recommendations are only indicative.

Platforms tested	FIA		
Antigens tested	Synthetic Brain Natriuretic Peptide-32 (BNP), Medix Biochemica 129-10 and synthetic BNP antigen, Medix Biochemica LA788.		
Product stability	temperatures affect the antigen bind	ife given on the first page is based on	
Miscellaneous	-		
References	-		



**Oy Medix Biochemica Ab** Klovinpellontie 3, Fl-02180 Espoo, Finland

medix@medixbiochemica.com www.medixbiochemica.com VAT reg. no. FI14631532

#### Legal disclaimer

Medix Biochemica products meet their specifications if transported, stored, and used according to the instructions. Medix Biochemica's products may not be used or reproduced without Medix Biochemica's written permission.