## Medix Biochemica

### Anti-h FDP LL01001 SPTN-5

Product overview				
Catalog number	C-40-0023			
Specificity	Antibody recognizes human fibrin degradation protein (FDP)			
Description	Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components.			
Product buffer solution	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN $_3$ as a preservative			
Shelf life and storage	18 months from manufacturing at 2–8 °C			
Analyte description	Fibrinogen has a molecular weight of 340,000 and is composed of three pairs of peptide chains (A $\alpha$ , B $\beta$ , C $\gamma$ ) connected by disulfide bonds. It transforms into fibrin under the action of thrombin. Both can be degraded into fragments such as X, Y, D, E, and D dimer by plasmin. Fragments X, Y, D, E, etc. are collectively referred to as FDP. FDP is a sensitive indicator of systemic or local intravascular coagulation and fibrinolysis.			
Parameters tested on each lot				
Product appearance	Liquid, may turn slightly opaque during storage			
Product concentration	5.0 mg/ml (+/- 10%)			

	<b>3</b> ( <b>1 1 1 1</b>
Immunoreactivity	80–120% compared to the reference sample
IEF Profile	5.8 - 6.9

≥ 95 %

### Kinetic parameters

Purity

Association rate constant	To Be Determined (TBD)

Dissociation rate constant TBD

Affinity constant TBD

Determination method Determination antigen



# Medix Biochemica

#### **Cross-reactivities**

Not Determined (N/D)

Epitope

Not Determined (N/D)

Pair recommendations

		DETECTION				
		LL01001	LL01002	LL01003	LL01004	LL01005
	LL01001	-	-	-	+	-
URE	LL01002	+	-	+	+	+
	LL01003	-	+	-	-	-
CAP'	LL01004	-	-	-	-	-
	LL01005	-	-	-	-	-

Please note that pair recommendations are based on results obtained by our laboratory. Equally good results may be obtained using other pairs and therefore these recommendations are only indicative.

Antigens tested	-		
Product stability	TEMPERATURE, TIME -70 °C, 21 days -20 °C, 21 days +4 °C, 21 days +35 °C, 21 days +45 °C, 7 days	RESULT OK OK OK OK OK	
	Stability testing is performed in the product buffer to see whether different temperatures affect the antigen binding, charge or composition of the antibody. Please note that the shelf life given on the first page is based on real time stability testing at 2–8 °C in the product buffer.		
Miscellaneous	-		
References	-		

